

Washable Isolation Gown

- Universal (One Size)
- 75 washes at infection control temperatures
- Comfortable, quick dry stretch cuffs safely cover the wrist area
- Wrap-around design for enhanced coverage with tie fastenings
- Carbon yarns reduce static build up
- Breathable fabric for comfort
- Designed with fewer seams for durability
- Stain resistant properties
- Fluid repellent fabric



Description: Washable Isolation Gown (Non-reversible)

Fabric: 99% Polyester, 1% Electrostatic Dissipative Carbon Yarns

Fabric weight: 105gsm

Dimensions: Centre back length: 123cm, Sweep: 153cm, Cuff: 8.9cm

Colour: Yellow

Testing: Fluid Repellent Fabric tested to AAMI Level 2
FDA registered

Laundering: Wash at Thermal Disinfection Temperature*
Do not use bleach
Do not exceed a maximum temperature of 90°C
Tumble Dry on a medium heat
Remove load immediately once goods are dry
75 washes at recommended care instructions

**The washing process should have a disinfection cycle in which the temperature of the load is either maintained at 65°C for not less than ten minutes or 71°C for not less than three minutes when thermal disinfection is used – (Department of Health - Choice Framework for local Policy and Procedures 01-04 – Decontamination of linen for health and social care: Management and provision)*

Please Note: All measures and weights may vary by plus or minus 5% and are based on new, unwashed products



wash



reuse



repeat



reduce
waste

Product Code: GOWN/PPE

Manufactured by



www.mipuk.co.uk

PRECISION
DENTAL PRODUCTS

www.precisiondental.co.uk
020 8236 0606

Care House, Chester Road, Borehamwood, Herts WD6 1LP
Precision Dental Products is a trading division of HealthCare Direct Ltd